

PIXX 1717 is a flat-panel digital x-ray detector with largest image coverage of 17" by 17" for every application with features not found on any other image device. PIXX 1717 is equipped with not only CsI but also premium Gadox scintillator for cost effective product line and it allows high X-ray sensitivity and high image quality. Thanks to 2nd generation AED technology which makes continuous and stable image acquisition without any loss of x-ray signal. On-board storage and Mobile Viewer allow users to use DR detector as like CR imaging plate. The detector provides other features that you love.



PIXX 1717

Flat-panel X-ray Detector



Product Options

Wired Standard (HDMI type) Data Link Cable
Wireless High Speed Wi-Fi, approx. 6 hours of battery life

Features

High Image Quality Faster Image Acquisition

AED Plus

Stand-alone Use On-board Storage

Web Viewer for Mobile Phone One-piece Carbon Body OLED Display for Status High SNR and Low Noise at lower dose

With largest image coverage 17" x17", less than 3~5 second

Second generation of continuous & accurate AED provides faster and efficient image acquisition Without workstation, detector works with x-ray source in portable & mobile environment

Protect any data loss by saving up to 100 images on the detector

Instant Viewer from any personal mobile devices

Full carbon fiber housing allows light and durable detector OLED window helps user to check the status of detector

Specifications

Application General Radiography

Sensor a-Si TFT array Flat Panel Detector

Conversion Screen CsI/GOS

Image Size 17 × 17 inch(427.25 mm × 427.25 mm)

Sensor Pixel 3,072 × 3,072 (9.4 M pixels)

 Pixel Size
 140 µm

 Spatial Resolution
 3.6lp/mm

 Image Data
 14 bit / 16 bit

Trigger Mode AED(Automatic Exposure Detection), Manual Trigger Calibration Mode ACC(Auto Calibration Control), Manual Calibration

Data Transfer TimeLess than 1 sec.Capture Cycle Time3 ~ 5 sec. approx.Wired InterfaceGigabit Ethernet

Voltage AC 100 ~ 240V, 50/60Hz I DC18V 3.5A

Wireless Standard IEEE 802.11ac

Battery Type Lithium Ion Polymer Battery
Dimensions (W x L x T) 461 mm x 460 mm x 15 mm

Weight 2.6 Kg (5.7 lb), 3.0 Kg (6.6 lb, with battery)









