



Compact with accurate and safe infusion
CareFusion Series
HK-100 series Infusion pump




Intuitive interface with colorful screen, all-in-one information at a glance .



Separate plastic covers for better waterproof.



Caregivers can easily install or remove pumps thanks to the rotatable pole clamp.



Detachable pump body for easy cleaning.



Anti-free flow design:
The Anti-free-flow(AFF) is aimed at preventing inadvertent excess infusion by automatically closing the tube with the anti-free flow clip when the door is opened.



Alarm management: With audible and visual alarm, caregivers can easily locate the alarm and warning.



Long-lasting battery: the battery life is long-lasting enough to last up to 5 hours.



High sensitive sensor guarantees the safe infusion by assist in detecting the risk and danger of air bubble and occlusion.



Anti-bolus helps prevent over-dose upon release of occlusion.