

### Features

- 10"LCD touch screen, folding 60 degree
- High sensitivity,10 chips 1MHz digital pulse transducer
- LE MEDICAL® DSP technology
- Twins monitoring function option
- Display waveforms and digital simultaneously
- FHR upper and lower limit alarm
- 6 kinds of printing mode, print paper size 112 mm
- Built-in high resolution and long-life thermal printer
- Automatic diagnosis, reforming Fischer score
- Doctor event mark, built-in expert dictionary
- 1200 hours data and waveform storage and review
- Wire connect to central monitoring system



### Technique Specification

#### • FHR

Ultrasound frequency:1.0MHz  
 FHR range: 40-240bpm  
 Ultrasound intensity:<5mW/cm<sup>2</sup>  
 Resolution:1 bpm  
 Accuracy: ±2 bpm

#### • TOCO

Range:0-100%  
 Resolution:1%  
 Fetal Movement  
 Manual fetal movement mark  
 Power Supply  
 AC 90V-240V,50/60 Hz  
 Battery: 7.4V 4000Mah