

As simple as it is precise

CareFusion Series HK-400 series Syringe pump

- Compact and light weight design allow easy moving and transferring.
- Automatically identifies various syringe types and sizes from 5ml to 60ml .
- Numeric and graphic pressure indicator helps caregivers identify occlusion in advance to prevent potential risks.
- deliver the precise medication dosage to even your tiniest patients. Flow rates deliver with +/- 2% accuracy as low as 0.1 mL/hr.
- Long-lasting battery: the battery life is long-lasting enough to last up to 5 hours.



Numeric keys combines button design facilitates fast programming.



Caregivers can easily install or remove pumps thanks to the integrated handle and rotatable pole clamp.



Pumps stack and lock together and can be expanded from 2 to 6(only for HK-400III).



With audible and visual alarm, caregivers can easily locate the alarm and warning.



Turn off the key sound and indicating light automatically with night mode on.



Hidden handle design helps to save room.

