

Features

-▶ 7"LCD screen, folding 60 degree
-▶ High sensitivity,10 chips 1MHz digital pulse transducer
-▶ LE MEDICAL® DSP technology
-▶ Twins monitoring function option
-▶ Display waveforms and digital simultaneously
-▶ FHR upper and lower limit alarm
-▶ 6 kinds of printing mode, print paper size 112 mm
-▶ Built-in high resolution and long-life thermal printer
-▶ Automatic diagnosis, reforming Fischer score
-▶ Doctor event mark, built-in expert dictionary
-▶ 1200 hours data and waveform storage and review
-▶ Wire connect to central monitoring system



Technique Specification

....▶ FHR

Ultrasound frequency:1.0MHz
 FHR range: 40-240bpm
 Ultrasound intensity:<5mW/cm²
 Resolution:1 bpm
 Accuracy: ±2 bpm

....▶ TOCO

Range:0-100%
 Resolution:1%
 Fetal Movement
 Manual fetal movement mark
 Power Supply
 AC 90V-240V,50/60 Hz
 Battery: 7.4V 4000Mah