

Multi Parameter

# Patient Monitor

PM12M

medimages



**Standard Parameter:** ECG, SpO2, NIBP, RESP, TEMP, HR

**Optional Parameter:** Touch Screen, Printer, 3/12-lead ECG, SD card, VGA output 2/4-IBP, EtCO2, AG module, Masimo/ Nellcor SpO2, Wall mount, Trolley.

## FEATURES:

- 1) 12.1" Color display
- 2) ECG, Spo2, NIBP, Resp, Temp
- 3) Display brightness adjustable.
- 4) Real time S-T segment analysis and ARR analysis.
- 5) SpO2 pitch tone variation and Drug Dose Calculation.
- 6) Parameter waveform and characters color selectable.
- 7) Multi-display: Standard, Large font, Trend Coaxis, Oxy CRG dynamic.
- 8) Standby mode for power saving and suspend monitoring.
- 9) 7-lead ECG waveforms display in phase.
- 10) Capture dynamic waveforms.
- 11) Built-in rechargeable lithium battery and battery volume checking.
- 12) Efficient resistance to interference of defibrillator and HF knife.

