



Our Ventilators are designed to provide critical respiratory support for patients in need. With advanced technology and user-friendly interfaces, our ventilators ensure precise and effective ventilation, enhancing patient outcomes and streamlining care.

15 inches colour touch screen

PATIENT TYPE:

Adult, Child, Infant

PARAMETER SETTING RANGE

Tidal Volume	: 2ml ~ 4000ml
Frequency	: 1 ~ 100 times/minute
Inspiration & Expiration	
Ratio	: 10:1 ~ 1:10
P _{supp}	: 0-80 cm H ₂ O
P _{insp}	: 5-80 cm H ₂ O
PEEP	: Off, 1-50 cm H ₂ O
Oxygen Concentration	: 21-100 vol%

INVASIVE VENTILATION MODES

VCV, PCV, VSIMV, PSIMV, CPAP/PSV, V+SIMV, BPAP, ASV, APRV, VS, PRVC

NON-INVASIVE VENTILATION MODES

PCV, SIMV-P, CPAP/PSV, APRV, BPAP, nCPAP, S/T

MONITORING PARAMETER

Tidal Volume	: 0-6000 mL
Frequency	: 0~200 times/minute
Inspiration & Expiration	
Ratio	: 150:1 ~ 1:150
Oxygen Concentration	: 18% ~ 100%
Leak	: 0% ~ 100%
PEEP	: ~120 cm H ₂ O
Minute Ventilation	: 0-100 L/Min

