



TS1 Transport Device
(Option)

Standard Configuration

Optional Configuration



Touch Screen



Fanless Design



Auto Brightness



WIFI



Plug and Play



Thermal Recorder



Color Waveform



Hot Key

1000

NIBP
Measurement
Review

120H

Trend
Review

1000

Alarm
Review

48H

Wave
Review

Rear View

- 1、Loudspeaker
- 2、AC Power Socket
- 3、Equipotential Grounding Terminal
- 4、Multi-function Interface
- 5、Network Interface
- 6、Dual USB interface
- 7、Dual USB interface
- 8、HDMI Interface

